## INSTRUCTIONS FOR USE

### ECG CABLES

**CAUTION**

**TEMPERATURE LIMITS**

**KEEP AWAY FROM DIRECT SUNLIGHT**

**PROTECT FROM MOISTURE**

**LOT CODE**

**QUANTITY**

**EXPIRY**

**DO NOT USE IF PACKAGING IS DAMAGED**

**MEDICAL DEVICE**

**NON-STERILE**

**MANUFACTURER**

**LATEX FREE**

Not made with natural rubber latex

**PHT FREE**

**READ INSTRUCTIONS FOR USE**

**REFERENCE**

**CONFORM TO MDR 2017/745**

ISO 13485:2016  
HS.NisshaMedical.com/IFU

## PATIENT ECG CABLES

### Please read carefully the following instructions

Failure to observe these precautions for use may lead to undesirable medical consequences.

### Important reminder:

This insert is designed to provide guidance for the use and handling of ECG cable which is designed for the measurement of the electrocardiogram in order to diagnose or monitor.

No reference is made to a specific medical technique. The manufacturer declines any responsibility for problems resulting from improper use of the product.

The classification rules may vary from country to country. According to EU MDR 2017/745 or Australian D Regulations, ECG Cables are products of Class I. According to the rules of FDA (USA) or Canada, ECG cables are **Class II**.

## I. IDENTIFICATION / FIELD OF APPLICATION

### IDENTIFICATION:

The **Nissha Medical Technologies ECG cable** is a crucial part of the system used to monitor or diagnose the cardiac electrical signal (ECG). It allows a safe connection between the electrodes on the patient and the electromedical device.

It is made of:

- a **trunk cable** fitted with two specific connectors:
  - o one of them connects the electromedical device
  - o and the other connects the leadwires, on the patient's side,
- three (or more) fixed (one-piece ECG cable) or detachable **leadwires**

The leadwires are fitted with means of connection (banana or pressure plugs or clips) allowing them to be connected to all of the various types of electrodes used.

### FIELD OF APPLICATION:

**Nissha Medical Technologies ECG cables** can be used on any kind of patient (adult, infant, paediatric and neonate) for both long and medium term ECG monitoring, resting ECG or for stress tests, on an outpatient basis, in the operating rooms, in Emergency sectors or in any other discipline requiring the measuring of an ECG.

There is a specific cable for each specific application site:

- The radiology room is for the use of X-Ray transparent leadwires
- The operating rooms are for the use of leadwires including HF protection

**Nissha Medical Technologies** offers a full range of **ECG cables** in its catalogue which can be viewed and downloaded on the company's website at: [www.NisshaMedical.com](http://www.NisshaMedical.com).

If you wish to place any subsequent orders, please use the code number shown either on the ECG cable or on its packaging.

If you require any further information about this product, please contact Nissha Medical Technologies or visit the company's website at [www.NisshaMedical.com](http://www.NisshaMedical.com).

## II. STORAGE / PACKAGING

### STORAGE

The storage conditions for Nissha Medical Technologies ECG cables are as follows:

- Ambient temperature: -40 to +70° C (-40° to +158° F)

### PACKAGING

Depending on the type (one-piece or detachable), **Nissha Medical Technologies ECG cables** are packaged either singly or in separate units (trunk cable, leadwires).

When it is not in use, the **ECG cable** must be stored in its original packaging in order to prevent any damage which might reduce its life time, performances and/or safety level.

## III. PERFORMANCES / RELIABILITY / SAFETY / COMPATIBILITY / MECHANICAL INTEGRITY / SYMBOLS / ALLERGICITY

### PERFORMANCES / RELIABILITY

Nissha Medical Technologies ECG cables are inspected both during and at the end of the manufacturing process according to technical protocols drawn up in line with current standards and regulations of this type of product.

The final results of the tests carried out on a group of representative ECG cables have been confirmed by an accredited Laboratory.

(LNE technical report no. 31273161.)

They have also undergone clinical testing and assessments.

(EC NO. DPL 034 PF 0018)

### SAFETY

**Nissha Medical Technologies ECG cables** are designed and manufactured in accordance with the general and special specifications of the relevant current European and/or international standards:

International standards IEC 60601-1 & 60601-2-25 / 60601-2-27

American standard AAMI EC 53

**Nissha Medical Technologies ECG cables** are part of the "applied part" to the patient as defined by the **EN IEC 60601-1** International safety standard.

The **ECG cables' safety class, type of protection (BF, CF) and degree of protection** against electrical shocks are closely linked to those of the electromedical device to which it may be connected to.

- **Please always remember to read the instruction manual for the electromedical device before using the ECG cable.**
- **Please keep the ECG cable at a distance as far as possible from any sources of electromagnetic radiation.**
- Please in the operating room make sure that the ECG cable (wholly or partly) is outside of the operating area in order to reduce the risk of burning the patient while the electrosurgical unit is in operation.
- In the operating rooms only use the specific high frequency protection ECG cable.

**Nissha Medical Technologies ECG cables** are designed to withstand repeated defibrillation shocks.

### PLEASE NOTE:

**Nissha Medical Technologies may not be held liable for any incidents which might occur to the patient or to medical staff caused by the presence of dangerous electrical currents coming from the electromedical device.**

### COMPATIBILITY

On the company's website ([www.NisshaMedical.com](http://www.NisshaMedical.com)), Nissha Medical Technologies offers its customers a downloadable document including information on the device's compatibility along with technical details about it.

### MECHANICAL AND ELECTRICAL INTEGRITY

The **Nissha Medical Technologies ECG cables** (conductors, connectors, cable, etc.) offer good mechanical resistance to traction and to flexion/bending in order to reduce the risk of damage during use, Nissha Medical Technologies has used materials offering high quality and excellent reliability. The moulded or assembled connectors are fitted with soft inserts designed to minimise the risk of the cable breaking at this point. The connector contacts are treated to reduce contact resistance.


(Test report no. CPB1097 – CPB1897/CPB2897)

The service life of an ECG cable depends on a large number of parameters. e.g.

- How many times it has been used
- At all times be sure to follow the safety instructions which you will find in the “PREVENTIVE MAINTENANCE” section below
- Being kept in a good state of hygiene

A regular visual and electrical inspection of the **ECG cable** will determine whether or not it needs to be replaced. (See preventive maintenance

**SAFETY SYMBOLS**

Shows that the ECG cables and accessories are fitted with special protection against electric shocks (including admissible leakage currents) and defibrillation shocks.

**The class and type of protection (BF, CF) against electrical shocks are defined by the type of electromedical device to which the Nissha Medical Technologies ECG cable is connected to.**

**ALLERGICITY**


The materials which come into contact with the patient used to manufacture Nissha Medical Technologies ECG cables have been subjected to allergicity tests. These tests have not shown the presence of any products which may lead to any intolerable allergic reaction.

**IV. INSTALLATION / USE / MAINTENANCE / HYGIENE / STERILISATION**  
**INSTALLATION**

Please follow these instructions for optimum installation and operation of the Nissha Medical Technologies ECG cable:  
(See also the instruction manual for the particular electromedical device and accessories)


1. If the leadwires are detachable, first connect the leadwires to the appropriate trunk cable corresponding to the correct location/colour code.

Example of detachable leadwires to be connected to the trunk cables:

Example of fixed (one-piece leadwires):

10 channels


Example of plugs fitted or moulded onto the leadwires:


Banana




Snap

Clip


2. Connect the plug of the **Nissha Medical Technologies ECG cable** to the appropriate electromedical device.  
 Plug Examples:


3. Depending upon the type of application, place the chosen electrodes onto the patient at the appropriate places for the measurement.

Code US

US COLOUR CODE (AHA)	
RA/ RIGHT ARM: <b>WHITE</b>	LA/ LEFT ARM: <b>BLACK</b>
LL / LEFT LEG: <b>RED</b>	RL / RIGHT LEG: <b>GREEN</b>
V1: 4th intercostal space, right edge of the sternum: <b>RED</b>	
V2: 4th intercostal space left edge of the sternum: <b>YELLOW</b>	
V3: on the 5th rib between C2 and C4: <b>GREEN</b>	
V4: 5th intercostal space, left medioclavicular line: <b>BLUE</b>	
V5: between C4 and C6 on the front left axillary line: <b>ORANGE</b>	
V6: on the left mid axillary line at C4 level: <b>VIOLET</b>	

Code EU

EU COLOUR CODE (ISO)	
RA/ RIGHT ARM: <b>RED</b>	LA LEFT ARM: <b>YELLOW</b>
LL/ LEFT LEG: <b>GREEN</b>	RL/ RIGHT LEG: <b>BLACK</b>
C1: 4th intercostal space, right edge of the sternum: <b>RED</b>	
C2: 4th intercostal space left edge of the sternum: <b>YELLOW</b>	
C3: on the 5th rib between C2 and C4: <b>GREEN</b>	
C4: 5th intercostal space, left medioclavicular line: <b>BROWN</b>	
C5: between C4 and C6 on the front left axillary line: <b>BLACK</b>	
C6: on the left mid axillary line at C4 level: <b>VIOLET</b>	

4. Connect the leadwires to the electrodes.
5. Fasten the **Nissha Medical Technologies ECG cable** nearest to the patient using the appropriate fastening system located on the lower side of the cable.
6. Start the electromedical device (see the instruction manual for the device).

**HOW TO USE:****GENERAL CONDITIONS:**

The conditions for using the Nissha Medical Technologies ECG cable are as follows:

- Ambient temperature: +10 to +40°C (+50 to +104°F)

**SPECIAL CONDITIONS:**

- **Do not use a cable or part of a cable if there is any risk to the patient (e.g. damaged insulating material).**

**PREVENTIVE MAINTENANCE:**

- If you need to take a cable out of the department for maintenance or verification, this becomes the responsibility of the user department to clean and disinfect the product before shipping or transportation.
- Make sure that the insulating materials (protective sheath on the connecting cable, interface cable, connectors) have not suffered any mechanical damage.
- Check the electrical continuity of the conductors monthly using an Nissha Medical Technologies physiological signal simulator. (Contact NISSHA MEDICAL TECHNOLOGIES to find out which is the most suitable tester for your needs)

**CORRECTIVE MAINTENANCE:**

- **There is no corrective maintenance for this product. If it is damaged it is to be discarded immediately.**

**HYGIENE:**

**Clean the ECG cable and its leadwires regularly using the method described below.**

**PLEASE NOTE:**

Do not allow any liquid to get inside the various connectors on the interface cable.

When using the Nissha Medical Technologies ECG cable in departments where a level of sterility must be adhered to, you must proceed to use the cleaning and/or disinfection methods described below.

**RE-PROCESSING INSTRUCTIONS OF ECG Cables (Acc. to EN ISO 17664:2004)**

<b>Re-processing limits</b>	NISSHA MEDICAL TECHNOLOGIES ECG cables are designed to support repetitive defibrillation shocks and cleaning and disinfection cycles. Any misuse of the cable may shorten its life time
<b>In the place of use :</b>	Remove any stains or dirt with a soft disposable cloth
<b>Isolation and transportation :</b>	Remember allways to reprocess the device after use.
<b>Cleaning preparation :</b>	No specific requirement.
<b>Automatic cleaning :</b>	Is strictly not applicable
<b>Manual cleaning :</b>	Clean the NISSHA MEDICAL TECHNOLOGIES ECG cable and its connecting cable (including the connector) with a cloth soaked in soapy water

<b>Disinfection:</b>	<p><b>Method A</b> (Taken from the study by the ANIOS Laboratory, no. 6416.94/0387)</p> <ul style="list-style-type: none"> <li>- fill a dip tank with a 0.5% HEXANIOS G+R solution</li> <li>- partially immerse the cable keeping the ends dry in order to avoid any electrical problems on the connectors</li> <li>- allow a contact time of 15 minutes</li> <li>- rinse the cable</li> <li>- dry the cable with absorbent paper</li> <li>- during the immersion rub the ends of the cable using a LINGET ANIOS wipe. (Replace the solution every 48 hours)</li> </ul> <p><b>Method B</b> (new cold method) (Taken from the study by the ANIOS Laboratory, no. 14496.02/052)</p> <ul style="list-style-type: none"> <li>- fill a dip tank with an ANIOXYDE 1000 activated preparation</li> <li>- but do not use before 30 minutes</li> <li>- partially immerse the cable protecting the ends in order to avoid any electrical problems on the connectors</li> <li>- allow a contact time of 10 to 30 minutes depending upon the required level of disinfection</li> <li>- regularly check the peracetic acid content using a test strip</li> <li>- and afterwards rinse with ANIOS mains water (pH = 7.3; TH = 48°f)</li> </ul>
<b>Drying</b>	Dry the cable with absorbent paper
<b>Maintenance, control and tests :</b>	Please check visually that the ECG cable is in working order
<b>Conditioning</b>	Please follow the standard protocol of packaging of your institution before sterilisation.
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>- Sterilisation time: 22 hours which takes into account 20 hours of gas exposure.</li> <li>- Sterilising agent: Carbon Dioxide + Ethylene Oxide (80/20%).</li> <li>- initial vacuum: -70 kPa.</li> <li>- relative humidity: &gt;60%.</li> <li>- temperature: 50 °C.</li> <li>- Ethylene oxide concentration: 530 g/m3.</li> <li>- pressure: 120 kPa.</li> <li>- final vacuum: -70 kPa.</li> <li>- aeration.</li> </ul> <p>(Test report 001: sterilisation test / ECG cables – RE/IP/BPF – Revision 0 dated 26/02/98)</p>
<b>Please Note :</b>	Never sterilise the cables and accessories in an “autoclave” using boiling water or steam.
<b>Preservation :</b>	Please read the conditions of storage in this document
<b>Complementary information :</b>	The above methods have been validated and are strongly advisable. Particularly we advise you strongly to use the recommended or equivalent products above and not to exceed the indicated time of application or immersion, otherwise the life time of the cuff might be reduced.
<b>Contact the manufacturer :</b>	See below or visit the website <a href="http://www.NisshaMedical.com">www.NisshaMedical.com</a>

**IV. GUARANTEE/LIABILITY**

Any unused NISSHA MEDICAL TECHNOLOGIES ECG cable kept in its original packaging and which has not suffered any visible damage is guaranteed by NISSHA MEDICAL TECHNOLOGIES for one year against hidden defects. NISSHA MEDICAL TECHNOLOGIES guarantees that the device complies with the specifications of the current safety and performance standards applicable to the ECG cables.

**PLEASE NOTE:**

 The class and type of protection (BF, CF) against electric shocks are defined by the type of electromedical device to which the cable ECG is connected.

Always read the instruction manual for the particular device and any accessories to be used before placing the device into service. NISSHA MEDICAL TECHNOLOGIES may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

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