INSTRUCTIONS FOR USE

HF CONNECTION CABLE FOR ELECTROSURGICAL PLATES

 TEMPERATURE LIMITS	 CAUTION
 PROTECT FROM MOISTURE	 NON-STERILE
 LOT CODE	 MANUFACTURER
 QUANTITY	 LATEX FREE <small>Not made with natural rubber latex</small>
 EXPIRY	 READ INSTRUCTIONS FOR USE
 MEDICAL DEVICE	 REFERENCE
 SEPARATE COLLECTION AS WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT	 EC ACCORDING TO EU MDR 2017/745

The classification rules may vary from country to country. According to the **EU MDR 2017/745** or the **Australian D Regulations**, HF connection cables are products of **Class I**. According to the rules of **FDA (USA) or Canada**, HF connection cables are **Class II**.

- The safety class, type of protection and degree of protection are defined by the electro-surgical device.
- Always read the instruction manual for the particular device and any accessories used with it before placing the HF connection cable into service.
- **Nissha Medical Technologies SAS** may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in the respective instruction manuals.
- Precautions must be taken in case the connection cable is disconnected on the electro-surgical device, in such a way as to prevent the metal parts of the connector which have become accessible from coming into contact with other live or earthed metal parts.
- Do not use a **HF connection cable** if there is any risk to the patient (e.g. damaged insulating material).
- If you need to take a cable out of the department to another department for maintenance or expertise, it is the responsibility of the department which has made use of the cable to clean and disinfect it before transferring or moving it.
- Never sterilize the **HF connection cables** for single use electro-surgical electrodes in an "autoclave" using boiling water or steam.
- The class and type of protection against electric shocks are also linked to those of the electro-medical device (Electrosurgery unit) to which the sensor and its cable are connected.
- Always read the instruction manual for the particular device and any accessories used with it before placing the device into service.
- **Nissha Medical Technologies SAS** may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

[HS.NisshaMedical.com/IFU](https://www.hs-nisshamedical.com/IFU)

HF CONNECTION CABLE FOR ELECTROSURGICAL PLATES

Please read the following information carefully.

Failure to observe these precautions for use may lead to undesirable medical consequences.

Important reminder:

This insert is designed to provide guidance for the use and handling of HF connection cables that are the link between the high frequency electro-surgical device and the single use "pre-gelled" electrode placed in contact with the patient.

No reference is made to a specific medical technique. The manufacturer declines any responsibility for problems resulting from improper use of the product.

I. IDENTIFICATION / FIELD OF APPLICATION

IDENTIFICATION:

The **HF connection cable for "electrosurgical plates"** provides the connection between the single use "**pre-gelled" electrode** placed in contact with the patient and the high frequency electro-surgical device. It is part of the circuit applied to the patient. Depending on the choice of neutral electrode used, it allows the application of safe HF energy to neutral electrode 1-zone or 2-zones. Within the limits of its specific features, it contributes to reduce HF shunt currents, which can cause risks to the patient, the user and other people.

The **HF connection cable** is made up of a conductor which is protected electrically by specific insulation and is fitted at each end:

- On the device's side, with a plug which is specific to each piece of high frequency electro-surgical device being used.
- On the patient's side, with a plug which is a connection system suitable for the type of single use neutral electrode being used.

The **HF connection cable** is available in two lengths of 3 m or 5 m.

FIELD OF APPLICATION:

HF connection cables can be used in all types of operating room equipped with a high frequency generator.

Nissha Medical Technologies SAS offers a full range of **HF connection cables for "electrosurgical plates"** in its catalogue (COMM/DOCU001/015) which can be viewed and downloaded on the company's website at [HS.NisshaMedical.com/IFU](https://www.hs-nisshamedical.com/IFU).

If you wish to place any subsequent orders, please use the code number shown either on your **Nissha Medical Technologies SAS HF connection cable** or on its packaging.

For further information about this product, please either contact Nissha Medical Technologies SAS or visit the company's website at [HS.NisshaMedical.com/IFU](https://www.hs-nisshamedical.com/IFU).

II. STORAGE / PACKAGING

STORAGE:

The storage conditions for **HF connection cables** are as follows

- **Ambient temperature: +10 to +50 °C (+50 to +122 °F)**

PACKAGING:

Nissha Medical Technologies SAS HF connection cables are packaged individually.

When it is not in use, the **Nissha Medical Technologies SAS HF connection cables** must be stored in its original packaging in order to prevent any damage which might reduce its service life, performances and/or safety level.

III. PERFORMANCES / RELIABILITY / SAFETY / COMPATIBILITY / MECHANICAL and ELECTRICAL INTEGRITY / ALLERGICITY

PERFORMANCES / RELIABILITY:

In order to ensure that the **HF connection cable** offers good resistance to traction and flexion, and to avoid untimely high frequency leakage currents as far as possible, **Nissha Medical Technologies SAS** uses materials and a manufacturing techniques which are both of high quality.

The molded connectors are fitted with soft inserts designed to minimize the risk of the cable breaking. **Nissha Medical Technologies SAS HF connection cables** are inspected both during and at the end of the manufacturing process according to technical protocols drawn up in line with current standards and directives on this type of product.

The final results of the tests carried out on a group of representative **HF connecting cables** have been confirmed by an accredited Laboratory.

(LNE technical report no. 5351010.)

They have also undergone clinical testing and assessments

SAFETY:

HF connection cables are designed and manufactured in accordance with the general and special specifications of the relevant current European, National and/or international standards:

(National, European and International Standards NF EN IEC 60601-1/NF EN IEC 60601-2-2)

(American standard AAMI HF18)

Nissha Medical Technologies SAS HF connection cables are part of the “**applied part**” to the patient as defined by the **IEC 60601-1 International safety standard**.

The **HF connection cable's** safety class, type of **protection (BF, CF, etc.) and degree of protection** against electric shocks are closely linked to those of the electromedical device to which it is connected.

Please always read the instruction leaflet for the electromedical device before using the **HF connecting cable** for electrosurgical plates.

HF connection cables are fitted with retractable tube connectors thus offering complete safety. There are no metal parts which can be in contact with the user during handling.

HF connection cables are designed to withstand repeated defibrillation shocks, disinfection and sterilization cycles. They have no accessible metal parts on the patient's side.

High and low frequency leakage currents, measured in accordance with the recommendations of the current standards which are applicable to this product, are below the authorized values.

**WARNING:**

The safety class, type of protection and degree of protection are defined by the electrosurgical device. Always read the instruction manual for the particular device and any accessories used with it before placing the HF connection cable into service.

Nissha Medical Technologies SAS may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in the respective instruction manuals.

Precautions must be taken in case the connection cable is disconnected on the electrosurgical device, in such a way as to prevent the metal parts of the connector which have become accessible from coming into contact with other live or earthed metal parts.

COMPATIBILITY:

On the company's website (HS.NisshaMedical.com), **Nissha Medical Technologies SAS** offers its customers a downloadable document including information on the device's compatibility along with technical details about it.

MECHANICAL & ELECTRICAL INTEGRITY:

In order to ensure that the HF connection cables (conductors, connectors, cable, etc.) offer good mechanical resistance to traction and to flexion and to reduce the risk of damage during use, **Nissha Medical Technologies SAS** has used materials offering high quality and excellent reliability. The moulded or assembled connectors are fitted with soft inserts designed to minimise the risk of the cable breaking at this point.

The service life of the HF connection cable depends on a large number of parameters. e.g.

- How many times it is used
- Adherence to preventive maintenance
- Being kept in a good state of hygiene

A regular visual and electrical inspection of the HF connection cable will determine whether or not it needs to be replaced. (See preventive maintenance)

The class and type of protection (BF, CF) against electrical shocks are defined by the type of electromedical device to which the Nissha Medical Technologies SAS ECG cable is connected to.

ALLERGICITY:

The insulating materials used to manufacture **Nissha Medical Technologies SAS HF connection cables** have been subjected to allergicity tests. These tests have not shown the presence of any products which might develop an intolerable allergic reaction.




INSTALLATION:

Please follow the instructions below for an appropriate installation of the HF connection cable for "monopolar surgical instruments":


(See also the instruction manual for the particular electrosurgical device and any other appliances being used)

- First of all connect the HF connection cable to the high frequency current electrosurgical device.
- Place "the single use neutral electrode" onto the patient at the appropriate location for the chosen application.
- As far as possible avoid any direct contact between the HF connection cable and the patient

Example of plug types on the electrosurgical device' side:

Example of plug types on the electrosurgical device' side:

HOW TO USE:

GENERAL CONDITIONS:

The conditions for using the **Nissha Medical Technologies SAS HF connection cables** are as follows:

- Ambient temperature: +10 to +40°C (+50 to +104°F)

SPECIAL CONDITIONS:

- **Do not use a HF connection cable if there is any risk to the patient (e.g. damaged insulating material).**

PREVENTIVE MAINTENANCE:

- Make sure that the insulating materials (protective sheath on the connecting cable) have not suffered any mechanical damage.
- Check the electrical continuity of the conductors (at least once per month).
(Contact Nissha Medical Technologies SAS to find out about the tester best suited to your own needs)
(Proper safety when using the connecting cable requires systematic inspections to be carried out least once per month or more, depending upon the frequency of use, disinfection or sterilisation.)
- Visual inspection (appearance of the insulating materials, the connectors' contact pins, etc.)
- Verification of the mechanical integrity of the connectors
- Check on the electrical continuity of the conductors
- Verification that the markings have remained properly in place.

WARNING:

- **If you need to take a cable out of the department to another department for maintenance or expertise, it is the responsibility of the department which has made use of the cable to clean and disinfect it before transferring or moving it.**

CORRECTIVE MAINTENANCE:

- **There is no corrective maintenance for this product. If it is damaged it must be discarded immediately.**

When using the **Nissha Medical Technologies SAS HF connection cables** in departments where a level of sterility must be adhered to, you must proceed to use the cleaning and/or disinfection methods described below:

RE-PROCESSING INSTRUCTIONS OF HF connections CABLES (Acc. to EN ISO 17664:2004)

Manufactured by: Nissha Medical Technologies SAS	Method:	Symbol:
Device(s): HF connection cable		
CAUTION	When cleaning or disinfecting the cables the plugs at the ends of the cables may become immersed in a liquid. Take care to dry those ends very carefully in order to avoid any electrical problem.	
Re-processing limits	Nissha Medical Technologies SAS HF Connection Cables are designed to support repetitive defibrillation shocks and cleaning and disinfection cycles. Any misuse of the cable may shorten its life time	
INSTRUCTIONS		
In the place of use	Remove any stains or dirt, particularly of organic origine, before they dry by cleaning the cable with a disposable soft cloth soaked in a warm soapy water or non aggressive detergent.	
Isolation and transportation	Remember always to reprocess the devices as soon as possible after use and to avoid any contamination do this on the spot.	
Cleaning preparation	No specific requirement.	

Manual cleaning	<ul style="list-style-type: none"> • Clean the Nissha Medical Technologies SAS HF Connection Cables (including the connector) with a cloth soaked in soapy warm water
Manual Disinfection	<p>(Taken from the study by the ANIOS Laboratory, no. 6416.94/0387)</p> <ul style="list-style-type: none"> • Fill a dip tank with a 0.5% HEXANIOS G+R solution • Partially immerse the cable keeping the ends dry in order to avoid any electrical problems on the connectors • Allow a contact time of 15 minutes • Rinse the cable • Dry the cable with absorbent paper • During the immersion rub the ends of the cable using a LINGET ANIOS wipe. (Replace the solution every 48 hours)
Drying after manual disinfection	Dry the cable with absorbent paper
Automatic cleaning	Not allowed
Maintenance, control and tests	Please check visually that the cable is in working order
Conditionnement	Please follow the standard protocol of packaging of your institution before sterilization.
Sterilization	<p>Process A:</p> <ul style="list-style-type: none"> • Sterilization time: 22 hours including 20 hours of gas exposure. • Sterilizing agent: Carbon Dioxide + Ethylene Oxide (80/20%) • Initial vacuum: - 70 kPa • Relative humidity: > 60 % • Temperature: 50 °C • Ethylene oxide concentration: 530 g/m³ <ul style="list-style-type: none"> • Pressure: 120 kPa • Final vacuum: - 70 kPa • Rinsing
Preservation	Please read the conditions of storage in this document
Complementary information	The above methods have been validated and are strongly advisable. Particularly we would advise you to strongly use the recommended or equivalent products above and not to exceed the indicated time of application or immersion; otherwise the life time of the cables might be reduced.
Contact the manufacturer	See below or visit the web site HS.NisshaMedical.com

V. GUARANTEE / LIABILITY

Any unused HF connection cable kept in its original packaging and which has not suffered any visible damage is guaranteed by Nissha Medical Technologies SAS for one year.

Nissha Medical Technologies SAS guarantees the device's conformity to the specifications of the current safety and performance standards applicable to it.

PLEASE NOTE:

- **The class and type of protection against electric shocks are also linked to those of the electromedical device (Electrosurgery unit) to which the sensor and its cable are connected.**
- **Always read the instruction manual for the particular device and any accessories used with it before placing the device into service.**
- **Nissha Medical Technologies SAS may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.**

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 **SANS LATEX**
Not made with natural rubber latex

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