



## Performance Plus<sup>™</sup>

Electrodes for Stress and Holter Testing

### FEATURES



#### HELPS REDUCE ARTIFACTS

Aggressive adhesive keeps electrodes in place on patient - reducing motion artifact and repeat testing.



#### HIGH CHLORIDE WET GEL

High chloride wet gel hydrates the skin and produces fast and clear traces.



#### Ag/AgCL SENSOR

Silver/Silver Chloride sensor provides accurate trace results.



#### HYPOALLERGENIC

Hypoallergenic adhesive, not made with natural rubber latex. Free from common allergens.

In critical diagnostic situations, choose PerformancePlus<sup>™</sup> Electrodes.

### PACKAGING

#### RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

#### CASE

FDA / GS1 Compliant Labels



#### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

#### INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | [HS.NisshaMedical.com](http://HS.NisshaMedical.com)

Patient Monitoring & Medical Surgical Consumables



Not made with natural rubber latex

Wet Gel






# PerformancePlus™ Electrodes




Electrodes for Stress and Holter Testing

In critical diagnostic situations, use PerformancePlus™ electrodes by Nissha Medical Technologies to ensure quality readings, fewer artifacts and reduce repeat procedures.




## A10005 - 1.5"/38.1mm - TEARDROP FOAM ELECTRODE

	Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	6500013V	A10005	5	5	100	600		 24 Hours
	6500020V	A10005NB	5	5	n/a	600		
	6500014V	A10005-1-60	1	60	n/a	600		
	6500016V	A10005-4NB	4	4	n/a	600		
	6500481V	A10005-7NB	7	7	n/a	700		
	6500018V	A10005-60	5	60	n/a	600		

## A10006-5 - 1.875"/ 47.625mm - TEARDROP FOAM ELECTRODE

	Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	6500024V	A10006-5	5	5	100	600		 24 Hours




## A10006-T - 2.25"/ 57.15mm - ROUND TAPE ELECTRODE

	Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	6500027V	A10006-5T	5	5	100	600		 24 Hours
	6500001V	A10006-1-60TS	1	60	n/a	600		

## A10007 - 2" / 50.8mm - ROUND CLOTH ELECTRODE

	Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	6500030V	A10007-5	5	5	100	600		 24 Hours

## A10012 - 1.875" x 1.25" / 47.625mm x 31.75mm - RECT. FOAM ELECTRODE

	Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	6500037V	A10012-5S	5	5	100	600		 24 Hours
	6500036V	A10012-4-60S	4	60	n/a	600		
	6500038V	A10012-5S	5	5	100	600		
	6500040V	A10012-60S	5	60	n/a	600		