

SilveRest™ Tab Electrodes

Diagnostic Resting EKG Tab Electrode

FEATURES



MINIMAL CLEANUP

Conductive adhesive solid gel requires minimal cleanup.



COMFORTABLE AND FUNCTIONAL

Pliable substrate conforms to the body.



CONVENIENT PACKAGING

Patented traction liner prevents tabs from sliding off the liner and ensures packs do not brick together for ease of removal and application.



HYPOALLERGENIC

Tab electrodes are not made with natural rubber latex.



Vermed SilveRest™ Tab Electrodes for resting EKG provide excellent adhesion and quality traces.

PACKAGING

RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

CASE

FDA / GS1 Compliant Labels



HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | HS.NisshaMedical.com

Patient Monitoring & Medical Surgical Consumables



Solid Gel



SilveRest™ Tab Electrodes

Diagnostic Resting EKG Tab Electrode

Vermed SilveRest™ Tab Electrodes by Nissha Medical Technologies for resting EKG provide excellent adhesion and quality traces. All of our tab electrodes meet or exceed the AAMI standard for disposable pre-gelled EKG electrodes.

A10022 - .984" x .90" / 24.99mm x 22.86mm - Ag/AgCL



Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6500043V	A10022	10	100	1000	4000		24 Hours

A10201 - 1.34" x .90" / 34.036mm x 22.86mm - C/AgCL



Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6500082V	A10201	10	100	1000	4000		24 Hours
6500083V	A10201-500	10	100	500	3000		24 Hours

A10023 - 1.125" x 1.438" / 28.575mm x 36.525mm - Ag/AgCL



Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6500044V	A10023	10	100	1000	4000		24 Hours
6502902V	A10023-1000BX	10	100	n/a	1000		24 Hours
6502971V	A10023-10	10	10		4000		24 Hours

A10202-10 - 1.297" x .864" / 34mm x 23mm - C/AgCL



Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6503179V	A10202-10	10	10	500	4000		24 Hours

A10027 - 0.5" x 0.9 / 12.7mm x 22.86mm - AG/AGCL TinyTab



TinyTab electrodes are specially designed to fit on smaller pediatric patients

Product #	Ref #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6502930V	A10027	10	100	n/a	1000		24 Hours