



# TenderTrode™

Tender Care for Tender Patients

## FEATURES



### INFANT BREATHABLE ELECTRODES

Designed for patient comfort. Gentle adhesive and low chloride hydrogel provide a conforming and non-irritating electrode. Excellent for use in APNEA monitoring and all other pediatric applications.



### SOFT CLOTH ELECTRODES

Adhesive solid gel electrodes with soft cloth material are available in various sizes to suit your needs.



### FOAM PEDIATRIC ELECTRODES


Our foam pediatric electrodes are available in 1" diameter to allow for easy placement and features solid gel that is gentle on the skin.

For use in NICU, Diagnostics, Apnea and PICU.

## PACKAGING

### RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes



### CASE

FDA / GS1 Compliant Labels



### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

1.800.669.6905 | [HS.NisshaMedical.com](http://HS.NisshaMedical.com)

Patient Monitoring & Medical Surgical Consumables



Solid Gel



Wet Gel




# TenderTrode™

Tender Care for Tender Patients

Vermed TenderTrode™ electrodes from Nissha Medical Technologies represent a complete selection of electrodes for infant and pediatric use. Each electrode is carefully designed with gentle, hypoallergenic adhesives and soft conforming materials to ensure comfort and effectiveness on young patients.

## A10001 - 1.625" / 41.275mm - ROUND TAPE ELECTRODE



Product #	Ref #	Per Strip	Per Pouch	Per Case	Wet Gel	Duration of Use
6500004V	A10001	3	3	600		
6500003V	A10000-60	3	60	600		72 Hours

## A10040 - 0.875" / 22.225mm - SQUARE CLOTH ELECTRODE



Product #	Ref #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6500060V	A10040-60	3	60	600		
						48 Hours

## A10041 - 0.875" x 1.3125" / 22.225mm x 33.34mm - RECT. CLOTH ELECTRODE



Product #	Ref #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6500061V	A10041-60	3	60	600		
6502767V	A10041-5	5	5	600		48 Hours

## A10042 - 1.375" / 34.925mm - SQUARE CLOTH ELECTRODE



Product #	Ref #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6500064V	A10042-60	3	60	600		
6500062V	A10042-10	5	10	600		48 Hours

## A10027 - 0.5" x 0.9 / 12.7mm x 22.86mm - AG/AGCL TinyTab ELECTRODE



Product #	Ref #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6502930V	A10027	10	100	1000		
						24 Hours