



VersaTrode™ Solid Gel

Adult General Purpose Electrodes

FEATURES



SOLID GEL ELECTRODES

Conductive adhesive solid gel improves adhesion and requires minimal clean up.



A10083 VERSATRODE™ REPOSITIONABLE

MR Conditional certified carbon snap and unique solid gel allows for the electrode to be repositioned up to 5 times.



A10071 VERSATRODE™ LEADLOCKING

Eliminate the need to tape down lead wires or replace displaced electrodes, saving time and money. Leadlocking electrodes safely secure leadwires, preventing artifact associated with movement.

Vermed VersaTrode™ electrodes are the ideal choice for all adult cardiac monitoring and diagnostic applications.

PACKAGING

RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

CASE

FDA / GS1 Compliant Labels



HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | HS.NisshaMedical.com

Patient Monitoring & Medical Surgical Consumables

MR Conditional



~~LATEX~~

Not made with natural rubber latex

Solid Gel



VersaTrode™ Solid Gel

Adult General Purpose Electrodes

Vermed VersaTrode™ electrodes are the ideal choice for all adult cardiac monitoring and diagnostic applications. Solid gel electrodes improve adhesion and reduce clean-up.

A10005SG - 1.5" / 38.1mm - TEARDROP FOAM ELECTRODE



Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case
6500021V	A10005SG-60	5	60	n/a	600
6500484V	A10005SG-30	1	30	n/a	600
6500484V	A10005SG-5NB	5	5	n/a	600
6502947V	A10005SG-3NB	3	3	n/a	600




A10033 - 1.875" / 47.625mm - TEARDROP CLOTH ELECTRODE



Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case
6500057V	A10033-1-60	1	60	n/a	600
n/a	A10033-5NB	5	5	n/a	600



A10083 - 1.25" x 1.25" / 31.75MM x 31.75mm - SQUARE CLOTH ELECTRODE RADIOLUCENT & REPOSITIONABLE



Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case
6500075V	A10083-5	5	5	n/a	600



A10096SG - 1.5" x 2.25" / 38.1mm x 57.15mm - RECTANGULAR FOAM ELECTRODE



Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case
6502754V	A10096SG-30	1	30	n/a	750

